


Lagging Block - Aluminium



All correspondence to:
Temprecord International Limited
PO Box 58 430, Botany, Manukau 2163, NEW ZEALAND.
Phone: (09) 274 9825, Fax: (09) 273 4020
Email/Web: info@temprecord.com, www.temprecord.com

Temprecord Lagging Block
(Fabricated to AS3864-1997 Standards)
Temprecord Part # 751100



NOTES:

- 1 The I.D. of the central hole should be such as to allow for the sliding fit of the probe at the operating temperature.
- 2 The thermal lagging block should be mounted at least 30 mm from any solid object and in the return air pathway where air circulation is fan forced.

All correspondence to:
Temprecord International Limited
PO Box 58 430, Botany, Manukau 2163, NEW ZEALAND.
Phone: (09) 274 9825, Fax: (09) 273 4020
Email/Web: info@temprecord.com, www.temprecord.com